

Parcan

**3% stabilized solution of sodium hypochlorite
for canal irrigation**

Composition

3% sodium hypochlorite, sodium chloride, sodium carbonate, solute qs ad 100g.

Indications

Canal irrigation for root canal preparation.

Properties

The root canal preparation, essential procedure for the success of the obturation, needs the use of irrigant solutions combined with a chemomechanical reaming.

PARCAN, a 3% sodium hypochlorite stabilized solution meets the requirements of the expected effects of an irrigation solution, i.e.:

- chemical effects by dissolution of the organic matrix of the dentin, allowing to flare the canal easily.
- the bactericidal potency of sodium hypochlorite.

PARCAN is a stabilized and buffered sodium hypochlorite solution ready for use.

Direction for use

Withdraw 10 ml of **PARCAN** solution in the syringe. Then carry out the canal irrigation by using appropriate needle to perfectly fit the syringe to avoid the leakage.

Warnings and precaution for use

So as to eliminate any risks of injecting the product beyond the apex it is strongly recommended :

- to check the apical closure (the use of the product for immature teeth or teeth with root resorption is contra-indicated)
- Use rubber dam kit to avoid contact of Parcan solution with soft tissues inside patient's mouth.
- Store away from heat and direct sun light.

Presentation

250 ml / 500 ml

For professional use only in the practice of dentistry.

Manufactured and Marketed in India by :
SEPTODONT HEALTHCARE INDIA PVT. LTD.
C-1/2, MIDC, Taloja, Taluka - Panvel,
Dist. Raigad - 410 208, Maharashtra, INDIA
Customer Care No. : 022-2740 1970 / 71
E-mail : info@septodontindia.com



R104105222

MADE IN INDIA

Customer : Septodont Healthcare	Location : NA	Date : 30 - 06 - 2022
Product Code : R104105222		
Product Name : Septodont PARCAN _Leaflet		Pharma Code:
Dimension : 270 X 100mm	Folding Size : 68 X 95 ± 1mm	
Substrate : 60 gms Maplitho Paper		
Leaflet Type : Insert		
Colour Front	K	
Non Printing Clours :		
Note - This approval will be considered for the final printing. Please recheck for corrections indicated earlier, In this proof also.		
Prepared By	Checked By	Proof Checked
		Approved By

PRABHAT PRINTER

Mob. : 9867390119